

PulmOne Installation Report

To be filled by a certified PulmOne Service Specialist

PO Number: Installation Date: Installer Name:

I. Product Details

Product Name:

Device S/N:

Tablet S/N:

HHU S/N:

II. Customer Details

Name:

Address:

City: Phone:

Email:

Customer's Staff	Name + Credentials	Phone / Cell	Email
Primary contact person			
Physician			
Practice Manager			
Professional #1			
Professional #2			
Professional #3			
Professional #4			
Professional #5			

III. Received Items

No.	Item	Y / N	Comments
1	MiniBox shipping box is sealed	☐ Y ☐ N	
2	Shipping contents intact:	☐ Y ☐ N	
	MiniBox device	☐ Y ☐ N	
	Handheld unit (Spirometer)	☐ Y ☐ N	
	Tablet	☐ Y ☐ N	
	Calibration filter	☐ Y ☐ N	
	Documentation USB	☐ Y ☐ N	
	Power and comm. cables	☐ Y ☐ N	
	Mobile cart	☐ Y ☐ N	
	Syringe	☐ Y ☐ N	

DLCO regulator	☐ Y ☐ N	
Filter pack - Qty: ____	☐ Y ☐ N	
3 Items supplied by customer	☐ Y ☐ N	

IV. Acceptance Tests

To be demonstrated and performed by PulmOne's Service Specialist

No.	Item	Y / N	Comments
1	MiniBox system is fully assembled and power is ON	☐ Y ☐ N	
	Touch screen turns ON and is responsive to touch	☐ Y ☐ N	
2	DLCO regulators and tank are connected to the MiniBox	☐ Y ☐ N	
	Regulator type: _____	☐ Y ☐ N	
	Pressure on gauges is: _____	☐ Y ☐ N	
	Pressure targets are marked	☐ Y ☐ N	
3	Cart is fully assembled	☐ Y ☐ N	
4	Supplied accessories assembled	☐ Y ☐ N	
5	Successful FVC calibration check (Error ____ %)	☐ Y ☐ N	
6	Successful DLCO calibration check	☐ Y ☐ N	
7	Successful FVC (performed by Service Specialist)	☐ Y ☐ N	
8	Successful SVC (performed by Service Specialist)	☐ Y ☐ N	
9	Successful LVM (performed by Service Specialist)	☐ Y ☐ N	
10	Successful DLCO (performed by Service Specialist)	☐ Y ☐ N	
11	Reports printed to a local printer	☐ Y ☐ N	
12	Printer brand & model: _____	☐ Y ☐ N	
13	All reports were saved to a USB flash drive	☐ Y ☐ N	
14	Additional setup information	☐ Y ☐ N	
15	Network drive connected successfully	☐ Y ☐ N	
16	Remote connection via TeamViewer was successful	☐ Y ☐ N	

V. Training

To be performed and confirmed by customer's staff

No.	Item	Y / N	Comments
1	Training conducted by PulmOne representative	☐ Y ☐ N	
2	Setup and configuration	☐ Y ☐ N	
3	Clinical settings	☐ Y ☐ N	
4	System settings	☐ Y ☐ N	
5	Calibration checks	☐ Y ☐ N	
6	Patient maneuvers demonstrated by PulmOne	☐ Y ☐ N	

7	Tests performed by staff on their patients:	☐ Y ☐ N	
	<i>FVC</i>	☐ Y ☐ N	
	<i>SVC</i>	☐ Y ☐ N	
	<i>MVV</i>	☐ Y ☐ N	
	<i>LVM</i>	☐ Y ☐ N	
	<i>DLCO</i>	☐ Y ☐ N	
	<i>Provocation</i>	☐ Y ☐ N	
8	Staff comfortable with all functions; enough training time given	☐ Y ☐ N	

VI. Post-Installation Activities

No.	Item	Y / N	Comments
1	Review service contact details	☐ Y ☐ N	
	<i>PulmOne website: www.pulm-one.com</i>	☐ Y ☐ N	
	<i>PulmOne toll-free service line: 1-844-PulmOne</i>	☐ Y ☐ N	
2	Create customer portal account for filter ordering (lead time 1-2 wks)	☐ Y ☐ N	
3	Review database backup option (customer's responsibility)	☐ Y ☐ N	
4	MiniBox Academy log-in	☐ Y ☐ N	

Customer Signature:

Date:

VII. Staff Checkup and Acceptance

Enter each staff member's name; mark Y or N for each item completed.

Item \ Staff member							
Volume calibration check	☐	☐	☐	☐	☐	☐	☐
DLCO calibration check	☐	☐	☐	☐	☐	☐	☐
FVC test	☐	☐	☐	☐	☐	☐	☐
SVC test	☐	☐	☐	☐	☐	☐	☐
MVV test	☐	☐	☐	☐	☐	☐	☐
LVM test	☐	☐	☐	☐	☐	☐	☐
DLCO test	☐	☐	☐	☐	☐	☐	☐
Provocation test	☐	☐	☐	☐	☐	☐	☐
Customer portal login access	☐	☐	☐	☐	☐	☐	☐
Filter ordering procedure	☐	☐	☐	☐	☐	☐	☐
MiniBox Academy signup	☐	☐	☐	☐	☐	☐	☐
Contacting Customer Service	☐	☐	☐	☐	☐	☐	☐
Setting up Settings	☐	☐	☐	☐	☐	☐	☐
Reports printing	☐	☐	☐	☐	☐	☐	☐
DLCO tanks replacement	☐	☐	☐	☐	☐	☐	☐
Signatures (initials)							

VIII. Comments / Outstanding Issues

PulmOne Service Specialist - Signature:

Date:

IX. Main Customer Person at Installation

Payment arrangements:

P.O.:

I confirm that I have been trained on the MiniBox(+).

I understand and am comfortable performing the following functions: general settings, clinical settings, report configuration, report printing / saving, calibration check, patient data entry, session start, and all patient test maneuvers; filter credits and updating sessions when a new box of filters is ordered.

I understand that the MiniBox+ has been validated to not require internal sanitization ONLY when the PulmOne viral-bacterial filter is used in conjunction with the MiniBox+. I also understand that the clinical accuracy of the MiniBox+ requires that a PulmOne viral-bacterial filter is used in conjunction with the MiniBox+. Use of the MiniBox+ without a PulmOne viral-bacterial filter may void your warranty.

I was also guided on how to access the customer portal and order mouthpieces and other supplies.

Customer Signature:

Date:

Print Name: