



MiniBox+

Cabinless Plethysmography



EASY

Easy to operate by anyone who performs spirometry. Fully automatic Lung Volumes and DLCO measurements. All results automatically generated.



FAST

Lung volumes with as little as 60 seconds of tidal breathing. Spirometry, lung volumes and DLCO in as little as 15-20 minutes.



ACCURATE

Clinically proven in trials with over 700 patients at 11 medical centers, based on comparative data from leading body plethysmography devices.

MiniBox+

Complete Pulmonary Function Testing (PFT)

Cabinless Plethysmography, Spirometry & Diffusion Testing

The MiniBox+ is the only portable, cabinless pulmonary function test (PFT) machine with gasless LVM plethysmography, as well as spirometry and diffusion testing (DLCO). With as little as 60 seconds of tidal breathing required for Lung Volume Measurements, this gasless, cabinless system provides maximum patient comfort.

FULLY COMPLIANT WITH ATS/ERS GUIDELINES

- ▶ **Spirometry**
Complete set of 30+ parameters including FEV1, FVC, SVC, MVV
- ▶ **Lung volumes - in as little as 60 sec**
Same volumes as body plethysmography including TLC, FRC, RV - with proven equivalency to body box¹
- ▶ **Diffusing capacity**
DLCO, VA, DLCO/VA - Single breath method



DLCO reference values may not be available for children under the age of 7

INNOVATIVE

Fully-automated Quality Control algorithms enable ease-of use and reliability


FULLY AUTOMATIC LVM AND DLCO

Manual measurement optional

ACCURATE & REPEATABLE

Automatically generated results

Z-SCORE



FULL SELECTION OF PREDICTED EQUATIONS

GLI reference values are standard for both Spirometry and Diffusion

PEDIATRIC USE

Measurements can be taken on children from age 5

SHORT TESTING TIME

Complete PFTs in as little as 15-20 minutes

GASLESS & CABINLESS PLETHYSMOGRAPHY

No sensation of claustrophobia; Can test obese, wheelchair-bound, or bedridden patients

MINIMAL TURNAROUND TIME BETWEEN PATIENTS

5-minute wipe-down, self-calibration, mobility for uninterrupted patient flow during room ventilation

CONNECTIVITY (OPTIONAL)

HL7 compliant, easily configurable, integrates with all major EMR systems


¹Data published in: CHEST 2021; 159(6):2356-2365

EASE OF USE

- ✓ Automatic zeroing of analyzers before each measurement
- ✓ No stabilisation time for plethysmography
- ✓ Warm-up time for analyzer less than 10 min
- ✓ Volume Gas graph for diffusion measurement
- ✓ Fast analyzer with response time t_{0-90} in less than 85 ms
- ✓ Time-to-test: as little as 15-20 minutes with pre- and post-spirometry

FLOW VOLUME

INCENTIVES



Complete mobility with the dedicated MiniBox+ cart

Ideal for:

Private Medical Practice

- ▶ Easy to operate by anyone on your medical office staff
- ▶ Fast testing – complete PFTs in only 15-20 minutes
- ▶ Comfortable for all your patients

Hospital Clinics

- ▶ Smallest footprint of any plethysmography device
- ▶ Fastest Complete PFT testing
- ▶ Most comfortable for your patients

Pediatrics

- ▶ Children can be tested from age 5
- ▶ Cabinless – not scary or claustrophobic
- ▶ Comfortable and fast measurements - increases compliance

Visit
our website:



Technical
specs:



MiniBox+ lung volumes are equivalent to the body box in a multi-center study

In a groundbreaking multi-center study published in **CHEST** in Feb 2021, lung volumes measured with the MiniBox+ were shown to be equivalent to body plethysmography ("body box") in both healthy participants as well as those with lung disease.

 266 subjects	 5 leading medical centers	 Compared with 3 different Body Box devices	 All ages, respiratory disorders, genders, severities	Study results MiniBox+ lung volumes were more reproducible and correlated better with body plethysmography than the other leading office-based techniques.
---	--	---	---	---

Mean Normalized Standard Difference (NSD)

	Helium Dilution	Nitrogen washout	Computerized Tomography	MiniBox+
%	15.3	14.2	13.1	8.9
Published studies	6	1	4	1

Mean Normalized Standard Difference (NSD) for TLC data when compared with body plethysmography

All tests are eligible for reimbursement*:

*subject to specific health plan policy

Code	Test
94010 or 94375	Forced vital capacity or Flow-volume loop
94060	Pre vs. Post bronchodilator - breathing capacity test
94200	Maximum breathing capacity
94726	Total lung capacity - determination of lung volume/plethysmography
94729	Diffusing Capacity
94618	Pulmonary stress testing (6-Minute-Walk-Test)

Coding should be based on medical necessity and is not intended to maximize reimbursement.

Enhance your diagnostics. With the MiniBox+.

ABOUT US

PulmOne is a global medical technology company dedicated to providing innovative solutions for pulmonary function testing in any clinical setting. The company's flagship product, the MiniBox+, is the only portable, cabinless pulmonary function test (PFT) device with gasless LVM plethysmography, as well as spirometry and diffusion testing (DLCO). The MiniBox+ meets all American Thoracic Society (ATS) and European Respiratory Society (ERS) guidelines.

PulmOne USA Tel: 1-844-785-6663 Email us: usinfo@pulm-one.com Visit us: www.pulm-one.com

PulmOne USA, 444 Madison Ave. 6th Floor, New York, NY 10022

PulmOne
PFT Outside The Box